






Brain-Computer Interfaces Lab Session P300 ERP

benjamin.wittevrongel@kuleuven.be






The P300 is a positive deflection of the EEG in response to a rare, attended stimulus (“oddball”)




central-parietal distribution



The P300 response can be used for BCI target selection, as it is only elicited by the intended target



In this session, we will solve a murder
by analyzing EEG for the P300 response



facts

butler has been found dead

9 suspects

9 possible murder weapons

9 possible murder locations

witness identified the details using a P300 BCI




your task

analyse the EEG

find the correct murderer, weapon and location

The witness identified the murderer using the following interface



Analyzing the EEG for the P300 response has following steps:



cut the EEG into epochs

[-0.1, 1.0] sec window time-locked to the onset of the stimulation

baselining

average pre-onset activity = 0


averaging per option

i.e., average epochs having the same marker

plot the ERPs

solve the murder

Analyzing the EEG for the P300 response:



<https://kuleuven.box.com/v/LabSessionsBCI>



Brain-Computer Interfaces Lab Session P300 ERP

benjamin.wittevrongel@kuleuven.be